

Determination of Interferon-Pathway-Related Gene Induction during RNAi



Jie Kang, Peter Hahn, Frank Narz, Cornelia Schmidt, Kirsten Haussuehl, Dong Liang,* Subu Yerramilli,* Eric Lader,* Dirk Löffert, and Wolfgang Bielke

QIAGEN GmbH, Hilden, Germany *QIAGEN Sciences, Germantown, MD, USA

Introduction

Experimental approaches using short interfering RNA (siRNA) molecules to specifically silence gene expression have become more widely used in recent years. As for all techniques in molecular and cellular biology, the importance of protocol optimization and the use of appropriate controls cannot be overestimated.

While the induction of interferon-pathway-related stimulatory processes during RNAi is usually attributed to long, double-stranded nucleic acids, it can also be a concern during transfection of short siRNA molecules or in other RNAi-related experiments. A number of parameters, including siRNA chemistry, transfection reagents, cell culture conditions, and siRNA concentration have the potential to cause interferon-pathway-related responses. These responses could interfere with the specific biological effect caused by the siRNA-mediated knockdown of the gene of interest, causing misleading results.

Here we describe assays for quantification of mRNAs which are derived from induced transcription of well characterized interferon-pathway-related genes, using quantitative, real-time RT-PCR. The assays are specific for STAT1, IL6, IFN α 1, IFN β 1, IFIT1, IFITM1, OAS1, or OAS2. We performed the assays after cells were transfected with functional amounts of synthetic siRNAs or with polyinosinic-polycytidylic acid (poly(I)-poly(C)), which is a long, double-stranded RNA that induces the interferon pathway. Although poly(I)-poly(C) transfection caused dramatic induction of several interferon pathway genes in these cells, as indicated by the assays, such induction was not seen when siRNAs from QIAGEN were transfected.

We conclude that these assays are valuable tools for optimization of experimental conditions and for use as controls, ensuring successful RNAi experiments.

A set of QuantiTect® Primer Assays for analysis of the interferon response

- QuantiTect Primer Assays are bioinformatically validated gene-specific primer sets for real-time, quantitative RT-PCR using SYBR® Green based detection (Figure 1).
- QuantiTect Primer Assays were used for analysis of the interferon-pathway genes STAT1, IL6, IFN α 1, IFN β 1, IFIT1, IFITM1, OAS1, and OAS2.

Principle of Real-Time PCR Using SYBR Green Detection




Figure 1 Primers bind to the nucleic acid template during the annealing step of PCR. Primers extend during the extension step and SYBR Green 1 dye binds to all double-stranded DNA molecules, emitting a fluorescent signal of a defined wavelength on binding.

Phosphorylation of eIF2 α induced by poly(I)-poly(C) but not by siRNA

- Phosphorylation of eukaryotic translation initiation factor 2A (eIF2 α) is a key event that occurs during the interferon response, resulting in nonspecific downregulation of protein synthesis.
- The level of phosphorylated eIF2 α (eIF2 α -P) in MCF-7 cells was measured after transfection of interferon-response-inducing poly(I)-poly(C), CDC2 siRNA, and a nonsilencing siRNA that contains 3 mismatches to CDC2 mRNA (CDC2-mut siRNA).
- The level of phosphorylated eIF2 α increased dramatically after transfection of poly(I)-poly(C). In contrast, very low levels of phosphorylated eIF2 α were observed after siRNA transfection (Figure 2).




Figure 2 MCF-7 cells were transfected with 1 μ g poly(I)-poly(C), 1 μ g CDC2 siRNA, or 1 μ g CDC2-mut siRNA (siRNA that contains 3 mismatches to CDC2 mRNA). After 4 and 24 hours, cells were analyzed by western blotting using antibody specific for eIF2 α and eIF2 α -P. **A** The resulting western blot is shown. **B** Quantification of western blot results.

Measurement of induction of interferon genes

- QuantiTect Primer Assays were developed to measure the expression of STAT1, IL6, IFN α 1, IFN β 1, IFIT1, IFITM1, OAS1, and OAS2 after transfection of poly(I)-poly(C) or a nonsilencing control siRNA into HeLa or MCF-7 cells.
- Interferon-pathway-related genes were upregulated after transfection of poly(I)-poly(C). This demonstrates that these assays are a highly effective tool for control experiments investigating the interferon response (Figure 3A and B).
- None of the interferon-pathway-related genes were upregulated after transfection of highly pure siRNA from QIAGEN using HiPerFect Transfection Reagent. This demonstrates that QIAGEN siRNA and HiPerFect Reagent do not cause nonspecific RNAi effects via induction of the interferon response (Figure 3).
- Strong induction of the interferon response by poly(I)-poly(C) but not by siRNA was also shown at the protein level (Figure 3C). siRNA transfection did not induce STAT1 protein expression, even using high siRNA concentrations. In contrast, poly(I)-poly(C) strongly induced STAT1 expression, even using as little as 1 ng.





Figure 3 **A** HeLa or **B** MCF-7 cells were transfected with 1 μ g poly(I)-poly(C) (abbreviated pl-C) or nonsilencing siRNA using HiPerFect Transfection Reagent. After 2 (MCF-7 only), 8, and 24 hours, gene expression was measured by quantitative, real-time RT-PCR using QuantiTect Primer Assays and the QuantiTect SYBR Green PCR Kit. **C** HeLa cells were transfected with a range of concentrations of poly(I)-poly(C), nonsilencing siRNA, or CDK4 siRNA using HiPerFect Transfection Reagent. After 48 hours, cells were analyzed by western blotting using STAT1-specific antibody and tubulin antibody (as an internal control). Untransfected (U) and mock transfected (M) cells were also analyzed.

Measurement of the interferon response in primary cells

- HUVEC were transfected with poly(I)-poly(C) or HP (HiPerformance) siRNA and the expression of IFN β 1 was determined after 72 hours, a typical timepoint for phenotypic knockdown evaluation. Knockdown was analyzed by quantitative, real-time RT-PCR.
- Poly(I)-poly(C) transfection induced IFN β 1 expression. siRNA transfection did not induce IFN β 1 expression, even after 72 hours in sensitive primary cells. Comparison with untreated cells showed that the viability of siRNA-transfected cells was not impaired.

Expression of Interferon Genes in Transfected Primary Cells



Transfection	C _t Value
HP siRNA	34.62
	33.78
	32.99
Poly(I)-poly(C)	19.78
	20.02
	20.05

Figure 4 HUVEC were transfected with 1 μ g poly(I)-poly(C) or HP siRNA (nonsilencing) using HiPerFect Transfection Reagent. Three replicate experiments were performed. After 72 hours, the level of expression of IFN β 1 was measured by quantitative, real-time RT-PCR using the QuantiTect Primer Assay and the QuantiTect SYBR Green PCR Kit.

Conclusions

- Various reports have shown that interferon-related pathways can be induced after siRNA transfection, leading to nonspecific effects and misleading results.
- We have developed QuantiTect Primer Assays for quantitative, real-time RT-PCR of the interferon-pathway-related genes STAT1, IL6, IFN α 1, IFN β 1, IFIT1, IFITM1, OAS1, and OAS2. These assays can be used as controls to determine the interferon response in RNAi experiments.
- Transfection with equimolar amounts of siRNA from QIAGEN or poly(I)-poly(C) using HiPerFect Transfection Reagent showed that the interferon response was not induced significantly by QIAGEN siRNA and HiPerFect Reagent.
- siRNA and HiPerFect Reagent are part of the range of RNAi solutions from QIAGEN that allow efficient knockdown without unwanted nonspecific effects.

Trademarks: QIAGEN®, QuantiTect® (QIAGEN Group); SYBR® (Molecular Probes, Inc.). siRNA technology licensed to QIAGEN is covered by various patent applications, owned by the Massachusetts Institute of Technology, Cambridge, MA, USA and others.

Purchase of QIAGEN products for PCR containing HotStarTaq DNA Polymerase is accompanied by a limited license to use them in the polymerase chain reaction (PCR) process for research and development activities in conjunction with a thermal cycler whose use in the automated performance of the PCR process is covered by the up-front license fee, either by payment to Applied Biosystems or as purchased, i.e. an authorized thermal cycler. The PCR process is covered by the foreign counterparts of U.S. Patents Nos. 4,683,202 and 4,683,195 owned by F. Hoffmann-La Roche Ltd.

QuantiTect Primer Assays are optimized for use in the polymerase chain reaction (PCR) process covered by patents outside the U.S., owned by F. Hoffmann-La Roche Ltd. No license under these patents to use the PCR process is conveyed expressly or by implication to the purchaser by the purchase of this product. Where the PCR process is covered by patents, a license to use PCR for certain research and development activities accompanies the purchase of certain reagents from licensed suppliers such as QIAGEN when used in conjunction with an authorized thermal cycler, or is available from Applied Biosystems. Further information on purchasing licenses to practice the PCR process where the process is covered by patents may be obtained by contacting the Director of Licensing, Applied Biosystems, 850 Lincoln Centre Drive, Foster City, California 94404 or the Licensing Department, Roche Molecular Systems, Inc., 1145 Atlantic Avenue, Alameda, California 94501.

© 2005 QIAGEN, all rights reserved.