


# Oxidation stability of sweet almond oil



Determination of the oxidation stability of sweet almond oil using the Rancimat method.

## Results

	Sample	Temperature	Induction time
—	Sweet almond oil	120 °C	4.1 h

# Method description

## Sample

Sweet almond oil

## Sample preparation

None

## Instrument

892 Professional Rancimat



## Parameters

Temperature	120 °C
Air flow	20 L/h
Evaluation	Induction time
Evaluation sensitivity	1.0
Stop criteria	End points
Sample amount	3 g
Measuring solution	60 mL deionized water