


# Automated Turnkey Solubility and Permeability Assays



## SPECTRAMAX APPLICATION NOTE #14

### INTRODUCTION

Millipore Corporation and MDS Analytical Technologies have partnered to provide platforms for solubility and permeability screening. The objective is to create automation-compatible 96- and 384-based systems with pre-validated methods, devices and analysis that can be run manually or in conjunction with standard laboratory automation. The integration with analytical methods whose throughput, robustness, and sensitivity are compatible with the needs of ADME screening significantly improves overall throughput.

The Millipore MultiScreen Solubility filter plate is a high-throughput, 96- and 384-well system to classify or quantify the aqueous solubility of compounds contained in typical libraries. The high correlation with the standard shake flask method, along with speed and low compound usage, gives this method an advantage over other methods, including nephelometry. The plate incorporates a low-binding membrane and low-binding plate materials to ensure accurate results.

MultiScreen filter plates used for PAMPA (parallel artificial membrane permeation assay) and MultiScreen Permeability filter plates increase throughput using non-cell-based, automation-friendly, 96- and 384-well assays that predict passive intestinal absorption. Non-cell-based assays streamline the ADME testing process, since they can be performed in less than 24 hours.

MDS Analytical Technologies SpectraMax<sup>®</sup> microplate readers for ADME screening are designed for high throughput in 96- and 384-well formats and offer software protocols

designed for the analysis and data reduction of assays performed using the MultiScreen filter plates. Protocols are available for permeability, qualitative solubility and quantitative solubility. The protocols control the instrument, perform data reduction and offer user-selectable criteria to flag the estimated 20% of compounds that are not amenable to UV-Vis detection. The readers allow parallel measurements and provide automation-friendly operation, with minimal methods development. The readers are available with FDA 21 CFR Part 11 compliance tools, as well as software and hardware validation tools.

### SOLUTION INCLUDES

MDS Analytical Technologies Microplate Readers for ADME screening

- SpectraMax<sup>®</sup> 190
- SpectraMax<sup>®</sup> Plus<sup>384</sup>
- SpectraMax<sup>®</sup> M2/M2<sup>e</sup>
- SpectraMax<sup>®</sup> M5/M5<sup>e</sup>
- FlexStation<sup>®</sup> 3

MDS Analytical Technologies SoftMax<sup>®</sup> Pro protocols for ADME screening

- MScreen Solubility Screen (qualitative solubility)
- MScreen Solubility Quantity (quantitative solubility)
- MScreen PAMPA (permeability)

Millipore Multiscreen Filter Plates (available from Millipore, Inc.)

- MultiScreen PBM Permeability Assay Plates  
Part Number: MAPBMN310
- MultiScreen Solubility Filter Plate (non-sterile)  
Part Number: MSSL BPC 10

### SALES OFFICES

#### United States & Canada

MDS Analytical Technologies  
Tel. +1-800-635-5577  
Fax +1-408-747-3601

#### Australia

MDS Analytical Technologies Pty. Ltd.  
Tel. +61-3-9896-4700  
Fax +61-3-8640-0742

#### Brazil

Molecular Devices Brazil  
Tel. +55-11-3616-6607  
Fax +55-11-3616-6607

#### China

Molecular Devices Beijing  
Tel. +86-10-6410-8669  
Fax +86-10-6410-8601

#### Molecular Devices Shanghai

Tel. +86-21-6887-8820  
Fax +86-21-6887-8890

#### Germany

MDS Analytical Technologies GmbH  
Tel. +49-89/96-05-88-0  
Fax +49-89/9-62-02-34-5

#### Japan

Nihon Molecular Devices Osaka  
Tel. +81-6-6399-8211  
Fax +81-6-6399-8212

#### Nihon Molecular Devices Tokyo

Tel. +81-3-5282-5261  
Fax +81-3-5282-5262

#### South Korea

Molecular Devices Korea, LLC  
Tel. +82-2-3471-9531  
Fax +82-2-3471-9532

#### United Kingdom

MDS Analytical Technologies Ltd.  
Tel. +44-118-944-8000  
Fax +44-118-944-8001

[www.moleculardevices.com](http://www.moleculardevices.com)

FLEXSTATION, MOLECULAR DEVICES, SOFTMAX, and SPECTRAMAX are registered trademarks of MDS Analytical Technologies (US) Inc. All other trademarks are the property of their respective owners.

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES.

©2009 MDS Analytical Technologies (US) Inc.  
Printed in U.S.A. 4/09 EP #0120-1473A

