

A complete workflow from sample preparation to analysis using SureSelect target enrichment system for Ion Proton semiconductor sequencing

Christian Le Cocq¹, Kyeong Soo Jeong¹, Arjun Vadapalli¹, Joseph Ong², Elin Agne³, Filip Karlsson³, Ashutosh Ashutosh¹, Francisco Useche¹, Jayati Ghosh¹, Henrik Johansson³, Scott Happe², Douglas Roberts¹, and Holly Hogrefe⁴
 Agilent Technologies, Diagnostics and Genomics Group ¹Santa Clara, California 95051 ²Cedar Creek, Texas 78612 ³Uppsala, Sweden ⁴La Jolla, California 92037





Abstract


Whole exome or targeted sequencing for protein-coding regions has provided a cost effective way to identify common and rare polymorphisms that are associated with Mendelian disorders and complex diseases. With increased capacity of semiconductor sequencing, highly multiplexed samples can be studied in a single sequencing run. However, a complete workflow processing raw DNA samples to identify DNA variants in target regions is not easily accessible. Here we describe an analysis workflow to study multiplexed samples in semiconductor sequencing for several target sizes: 50Mb (Human All Exon), 3.2Mb (Human Kinome) and a 1Mb custom design. The workflow includes library preparation, SureSelect target enrichment, semiconductor sequencing, and variant calling with SureCall software (beta version). Improved and simplified steps for library preparation and target enrichment maximize multiplexing and produce consistent results in the Ion Proton sequencer. Sequencing output can be easily analyzed, visualized and summarized in a report with SureCall which is optimized for use with Agilent's target enrichment system. We demonstrate high capture efficiency, uniformity, and reproducibility of enrichment. The results from different capture sizes show comparable high performance regardless of various targeted regions. The combination of efficient target enrichment system, semiconductor sequencing, and SureCall software provides a fast and convenient tool to assess DNA variants in genomic regions of interest.

Materials & Methods

SureSelect Target Enrichment Workflow



SureCall Workflow



All Exon V5 Capture Performance

SureSelect with Ion Proton shows excellent capture and uniformity.




Figure 1. Summary of typical capture performance. SureSelect was performed with Hapmap DNA (NA12878) and Exon V5 baits. The vertical axis shows the percentage of on-target or coverage. Each sample includes 4Gb of sequencing. Coverage shows the percentage of targeted bases with triplicates.

Read Depth Distribution




Figure 2. Read depth distribution. The uniquely mapped reads show a uniform read depth distribution. The cumulative percent of reads is shown by an orange line.

SNP concordance (All Exon V5)

Evaluation sites	49113
Overlapping sites with db135	39608
Concordant sites	39435
Novel sites	9433
Concordant rates (%)	99.4

Table 1. SNP concordance with dbSNP135. The exon data show high concordance with previously reported SNPs. The amount of sequencing used for comparison is 4 Gb. SNP calling was performed with Genome Analysis Toolkit.

Kinome Capture Performance




Figure 3. Capture performance with SureSelect Kinome baits of 3.2 Mb. 8 Hapmap DNA (NA12878) samples were analyzed in one Ion Proton run with multiplexing. Each sample was normalized with 0.32 Gb of sequencing. Coverage represents the percentage of targeted bases with at least 1 (blue), 10 (orange) and 20 (green) reads.

Custom Capture Performance




Figure 4. Capture performance with SureSelect Custom baits (1 Mb). 16 Hapmap DNA (NA12878) samples were analyzed with multiplexing. The results are shown with 0.1 Gb of sequencing. The vertical axis shows the percentage of targeted bases with at least 1 (blue), 10 (orange) and 20 (green) reads.

SureCall results



Figure 5. Triage view in SureCall v1.1 highlighting dbSNP record for SNP rs6829402 on the left and indel rs33992431 on the right for NA12878. The Triage view allows users to view variants and the reads supporting the call in a single view along with annotations from public databases.



Figure 6. Corresponding dbSNP records shown for SNP (left) and indel (right).

Coefficient of Variation	Percentage of supporting reads	Total Number of Variants	Percentage of reported Indels
no filter	no filter	10751	87
3	75	4326	66
1	75	3718	61
0.5	75	1512	4

Table 2. Indel Filter performance in SureCall v1.1 on Ion Proton SureSelect Kinome data. The higher rate of indel calls on the Ion platform is primarily due to inaccurate flow calls over homopolymer regions. The number of indels reported in the Triage view is controlled by the parameters coefficient of variation and percentage of supporting reads. The coefficient of variation parameter is defined as the standard deviation over the mean indel size. The percentage of supporting reads is defined as the percentage of reads supporting the indel call.

Conclusions

- Agilent's SureSelect Target Enrichment for the Ion Proton Platform provides a comprehensive, efficient, robust, and cost-effective means to sequence subsets of the human genome.
- Different capture sizes show comparable high performance regardless of various targeted regions.
- High reproducibility of enrichment, depth distribution, and sequence coverage from multiplexed sequencing.
- Excellent concordance with known dbSNP.