


High-throughput analysis of protein formulations by DLS: thermal stability, colloidal stability and a thermal anomaly in the colloidal stability parameter D_1

Sophia Kenrick and Daniel Some, Wyatt Technology Corporation

1. Introduction

High-Throughput Screening by Dynamic Light Scattering (HTS-DLS) is a versatile tool for characterizing various aspects of protein stability *in parallel*:

- Thermal / conformational stability (T_M , T_{onset})
- Colloidal stability (D_1 or k_D)*
- Actual aggregation state



D_1 undergoes a dramatic, previously unobserved modulation around the thermal transitions.

Changes to in D_1 are observable several degrees earlier than T_{onset} . We speculate that this reflects an unfolding transition.


*we prefer 'D₁' to 'k_D' in order to avoid confusion with other 'K-D' meanings and to open a pathway to higher-order terms.

2. Model System

mAb1

- pI = 8.9 by electrophoretic light scattering (Möbius®, Wyatt Technology)
- R_H : 4.7-4.8 nm at 25°C, c→0
- pH conditions: 8.5 ($Z_{DHH}=4$), 9.5 ($Z_{DHH}=-2$)
- Concentrations: 0.47, 0.94, 1.88, 3.75, 7.5, 15 mg/mL
- 3 replicates per pH & concentration
- Wells capped with silicone oil
- Measured every 0.5°C over 25°C - 85°C

3. Analysis



T_M is the midpoint of a thermal transition in the hydrodynamic radius R_h (calculated from the mutual diffusion constant D_m under dilute conditions).

T_{onset} is the temperature at the onset of aggregation, determined from R_h .

D_1 is determined from the concentration dependence of D_m :

$$D_m = D_0(1 + D_1c + D_2c^2 + \dots)$$


D_1 is closely related to the second virial coefficient A_2 (B_{22}) and is an indicator of colloidal interactions.

HTS-DLS of T_M/T_{onset} , D_1 and aggregation size distributions under multiple formulation buffer conditions, can be accomplished simultaneously in situ in standard low-volume microtiter plates with a DynaPro® Plate Reader II.

All parameters determined from cumulants.

4. DynaPro® Plate Reader II


- Standard 96, 384 or 1536 microwell plates
- Sample stays in wells – no carryover
- Measurement time = 4-10 s/well, typical
- Temperature ramps 4-85°C
- Camera for troubleshooting bad data



5. Thermal Stability


- pH 8.5** • $T_{onset}(c)$:
- ~56°C (0.47 mg/mL)
 - ~53°C (15 mg/mL)

- Degree of aggregation depends on concentration
- Quasi-stable state between 60°C-70°C




- pH 9.5** • $T_{onset}(c)$:
- ~55°C (0.47 mg/mL)
 - ~52°C (15 mg/mL)

- Degree of aggregation depends on concentration
- Stable state above ~ 61°C



6. Aggregation States


- HTS-DLS indicates vastly different aggregation states, not expected to be identified by DSC or DSF signal loss
- Does form of $D_1(T)$ point to aggregation pathway?



8. Colloidal Parameter D_1


- pH 8.5** • ' $D_1(T)$ ' – transition midpoint at T_{onset} ;
- 'Stabilizes' with R_h
 - May not have much meaning

$T > T_{onset}$ (new value just a consequence of earlier concentration-dependent aggregation?)



- pH 9.5** • $D_1(T)$ – negative peak at T_{onset} !
- 'Stabilizes' with R_h , 2nd transition ~75°C

• After transition – true interaction between aggregates or consequence of earlier aggregation?



9. Conclusions & Future Studies

High-throughput analysis by the DynaPro DLS plate reader provides parallel analysis of *multiple* stability-indicating parameters for **screening of developability and optimal formulations**.

A novel indicator, $D_1(T)$, provides new insights into the **interplay of colloidal and conformational stability**. The form of $D_1(T)$ may point to aggregation mechanisms and rate-limiting parameters.

Better understanding of the significance of $D_1(T)$ will require further study including tests of reversibility.

The nature of $D_1(T)$ above the transition may be further elucidated by a dilution series performed on aggregated samples in the context of HTS-DLS.