


# LC/MS/MS quantitation of eicosapentaenoic acid in rat plasma

Rachel Sun, Jordan Nally, Tim Shoaf  
BASi, 2701 West Lafayette, IN 47906

## Introduction and Purpose

Eicosapentaenoic acid (EPA) is an omega-3 fatty acid, which has been documented to be important on biochemical and physiological processes including inflammatory response, blood clotting, and the immune system. Eicosapentaenoic acid (EPA) is different forms in the blood circulation, such as triglyceride, ethyl ester, or free acid. In this work, we presented an LC/MS/MS method to quantitate the total eicosapentaenoic acid in rat plasma

## Structure of EPA



## Extraction procedure

Transfer sample/Cal/QCs/blanks to 96 well plate  
Add internal standard working solution  
Add 200 µL of hydrolyzation solvent ACN/6N HCl (90:10)  
Hydrolyze the samples at 100 °C for 1 hour and 15 min  
Cool down the plate to room temperature. Centrifuge the plate briefly.  
Add 650 µL of hexane for liquid-liquid extraction  
Transfer the supernatant to the clean plate with 10 µL DMSO in each well.  
Evaporate to dryness and add reconstitution solutions  
Vortex and centrifuge.  
Seal the plate and inject samples.

## LC/MS Conditions


Column: Kinetex, XB-C18, 50\*2.1, 2.6 µm  
Mobile Phase A: 10% ACN/90% water,  
Mobile Phase B: 90% ACN/10% water, 10 mM Ammonium acetate.  
Run time: 3.5 minutes  
TSQ Quantum mass spectrometer  
Source: HESI  
Ion mode: Negative  
Resolution (Q1/Q3): Unit/Unit  
MRM masses (\*): 301/257 (EPA)  
306/262 (EPA-d5)

## Results

### Calibration statistics

Assay Date	Analytical Run Number	STD 0.500 0.500 µg/mL	STD 1.000 1.00 µg/mL	STD 2.000 2.00 µg/mL	STD 5.000 5.00 µg/mL	STD 20.0 20.0 µg/mL	STD 50.0 50.0 µg/mL	STD 100 100 µg/mL	STD 250 250 µg/mL
18-Mar-2013	1	0.526	1.08	2.18	5.21	21	53.4	104	277
18-Mar-2013	2	0.481	0.888	1.82	4.9	18.4	49.3	92.7	224
19-Mar-2013	3	0.584	1.11	2.13	5.17	19.7	50.2	101	253
19-Mar-2013	4	0.485	1.04	2.03	5.35	20.5	53.2	104	263
19-Mar-2013	5	0.519	0.985	1.77	4.98	18.8	47.8	99.9	233
19-Mar-2013	6	0.412	1.01	1.93	5.04	19.4	51.4	100	248
20-Mar-2013	1	0.598	0.971	1.98	5.09	19.5	51.2	105	243
20-Mar-2013	2	0.522	0.949	1.89	5.27	19.7	50.4	98.8	256
20-Mar-2013	3	0.51	0.899	1.74	5.01	20.1	51.2	105	237
20-Mar-2013	4	0.526	1.02	1.98	5.12	20.3	50.6	108	246
Mean		0.504	0.997	1.94	5.1	19.7	50.7	102	248
S.D.		0.0454	0.0708	0.133	0.131	0.714	1.65	3.98	13.9
%CV		9	7.1	6.9	2.6	3.6	3.3	3.9	5.6
%Bias		0.8	-0.3	-3	2	-1.5	1.4	2	-0.8
n		12	12	12	12	12	12	12	12

### Chromatograms




Proxy Blank-

LLOQ

ULOQ

### Typical calibration curve



### QC statistics

Run Date	Curve Number	QC 2.00 2.00 µg/mL	QC 80.0 80.0 µg/mL	QC 200 200 µg/mL
18-Mar-2013	1	2.02	83.4	221
18-Mar-2013	2	1.78	77.2	200
18-Mar-2013	3	2.07	80.9	206
19-Mar-2013	1	2.08	78.9	208
19-Mar-2013	2	2.11	77.2	206
19-Mar-2013	3	1.93	76.1	199
19-Mar-2013	4	2.02	77.4	199
19-Mar-2013	5	1.93	80.1	263
20-Mar-2013	1	2.17	82	221
20-Mar-2013	2	1.87	81.4	194
20-Mar-2013	3	2.16	82.2	185
Mean		2	79.7	208
S.D.		0.121	2.96	19.8
%CV		6.1	3	9.5
%Theoretical		100	98.6	104
%Bias		0	-0.4	4
n		12	12	12

## Discussion

proxy matrix for endogenous compound: parallelism test

### Calibrators in proxy matrix

Assay Date	Analytical Run Number	STD 0.500 0.500 µg/mL	STD 1.00 1.00 µg/mL	STD 2.00 2.00 µg/mL	STD 5.00 5.00 µg/mL	STD 20.0 20.0 µg/mL	STD 50.0 50.0 µg/mL	STD 100 100 µg/mL	STD 250 250 µg/mL
21-Feb-2012	1	0.547	1.02	1.95	5.21	19.7	49.3	111	239
		0.46	0.963	1.96	4.96	19.5	47.4	105	249
Mean		0.504	0.992	1.96	5.09	19.6	48.4	108	244
S.D.									
%CV									
%Bias		0.8	-0.8	-2	1.8	-2	-3.2	9	-2.4
n		2	2	2	2	2	2	2	2

### QCs in plasma and proxy matrix

Run Date	Curve Number	Proxy matrix			Plasma Matrix		
		QC 2.00 2.00 µg/mL	QC 80.0 80.0 µg/mL	QC 200 200 µg/mL	QC 2.00 PL 8.00 µg/mL	QC 80.0 PL 86.0 µg/mL	QC 200 PL 206 µg/mL
21-Feb-2012	1	2.08	76.4	199	9.25	85.8	199
		2.5	78.9	198	8.36	87	194
		2.09	75.2	190	8.24	85.7	195
		1.95	78.1	190	8.53	86.5	202
Intran Mean		2.16	77.2	191	8.57	86.3	195
Intran SD		0.239	1.67	1.73	0.472	0.614	7.87
Intran %CV		11.1	2.2	0.9	5.5	0.7	4
Intran %Bias		8	-3.6	-4.5	7.1	0.3	-5.3
n		4	4	4	4	4	4

### Hydrolysis of EPA conjugates

Majority of EPA is found in plasma phospholipids and lipoprotein complexes. Plasma samples are hydrolyzed at acidic condition to release EPA from conjugates. The hydrolysis was carried out in a 96 well plate with an in house made clamp to prevent acetonitrile evaporation.

## Conclusions

A fast, simple, specific, and robust method was developed with excellent linearity, precision, and accuracy for the quantification of EPA in rat plasma. The method has been used successfully to assay samples from non GLP studies.

## Contact Information

BASi  
2701 Kent Ave  
West Lafayette, IN 47906  
www.basinc.com

Contact:  
Rachel Sun  
rsun@basinc.com