


Specificity and functionality of microRNA inhibitors

Barbara Robertson, Andrew Dalby, Jon Karpilow, Anastasia Khvorova, Devin Leake and Annaleen Vermeulen Dharmacon, now part of GE Healthcare, 2650 Crescent Drive, Suite #100, Lafayette, CO 80026, USA



Introduction

Highly potent microRNA (miRNA) inhibitors are valuable tools for elucidating the roles of miRNAs and their targets. Although it is known that pairings between endogenous miRNAs and their natural targets in animals generally involve base pair mismatches, studies of how mismatches between endogenous miRNAs and artificial inhibitor targets might affect inhibitor specificity and functionality have been very limited. Using a Luciferase reporter system we have investigated the specificity of miRNA inhibitors. We first confirmed significant levels of cross-reactivity among the closely related let-7 family members. Subsequently, a systematic study of mismatches incorporated into inhibitors of single-family-member miRNAs identified two regions that are important for overall inhibitor functionality. Our findings indicate that features important for natural miRNA target recognition also appear to be important for inhibitor specificity. Understanding the specificity of inhibitors allows for better interpretation of inhibitor activity in endogenous systems.


Systematic testing of mismatched inhibitors





Mismatches affect expression of an endogenous target with similar positional dependence



Cross-reactivity is evident among let-7 microRNA family member inhibitors



Mismatched inhibitor target sites show that position of mismatch affects functionality





Examples of variation in sequence identity among microRNA family members and a non-family member with matching seeds

miRNA name	Mature sequence from 5' to 3' end
Family members with sequence identity in both the seed region and the 3' region	
hsa-miR-15a ¹	UAGCAGCACAUAAUGGUUUGUG
hsa-miR-15b ¹	UAGCAGCACAUCAUGGUUUAACA
Family members with sequence identity in the seed region but very limited identity in the 3' region	
hsa-miR-15a ¹	UAGCAGCACAUAAUGGUUUGUG
hsa-miR-16 ¹	UAGCAGCAGGUAUAUUGGGCC
Non-family members with sequence identity both in the seed region and in the 3' region	
hsa-miR-15a ¹	UAGCAGCACAUAAUGGUUUGUG
hsa-miR-497 ²	CAGCAGCACACUGUGGUUUGU

Notes: The seed region is in *italics*, the 3' region (position 13-18) is bold
¹Member of the miR-15 family
²Member of the miR-497 family

Sequence similarity between miR-21 and miR-122 does not cause cross-reactivity between inhibitors



Conclusion

- Mismatch position in the target region of the inhibitor affects inhibitor efficacy. The seed region (position 3-8) and a 3' region (positions 13-18) appear to be equally important for inhibitor functionality and specificity.
- Results of the mismatch study can help identify potential cross-reacting microRNAs.
- Future application of these results may help design more effective inhibitor experiments: For miRNAs whose inhibitors are expected to cross-react, use of a pool of inhibitors may give more consistent results. For miRNA family members whose inhibitors may not all be expected to cross-react, use of a pool of inhibitors may reveal a phenotype otherwise obscured by redundancy of miRNA function.

Methods

- All inhibitors are fully 2'-O-methylated molecules synthesized at at Dharmacon according to the miRIDIAN Hairpin inhibitor design.
- Dual-luciferase reporters were derived from the psiCHECK-2™ vector (Promega) that contains both firefly (FLUC) and Renilla (RLUC) luciferase genes. The desired microRNA target was cloned into the 3' UTR of the RLUC gene generating a cleavage assay where the mature miRNA is fully complementary to a single site. Results shown are averages from triplicate wells, normalized to corresponding controls, then expressed as fold-inhibition relative to negative control. Error bars are +/- one sample standard deviation of the original triplicate data, scaled for all subsequent calculations.
- Levels of aldolase A (ALDOA) messenger RNA (mRNA) (a known target for miR-122) were determined using the QuantiGene Assay (branched DNA; Affymetrix, Panomics products).

gelifesciences.com/dharmacon

Dharmacon™
part of GE Healthcare

