

Application of Solid Phase Micro Extraction (SPME) and On-Fiber Derivatization

to Detect Precursors of Biomolecules in the study of Proto-metabolism


— プロトメタボリズム研究における生体分子前駆体検出のためのマイクロ固相抽出・誘導体化法 —

○Y.Hongo¹, K.Yamada², H.Hübschmann³, R.Yi¹

(WPI-ELSI, Tokyo Institute of Technology¹, AMR, Inc.², CTC Analytics AG³)


Keywords: proto-metabolism, GC-MS, SPME, on-fiber derivatization

Prebiotic organic syntheses on the primitive early earth (4.5-2 Ga)



In the studies of origins of life, plausible synthetic reaction pathways to form biological relevant molecules, or building blocks of proto-bi-layer lipids, protein, RNA, and DNA in the primitive environment have been investigated intensively. We are focusing on studying syntheses of primitive molecular set, including **Sugar, amino acid, and their analogues**, which were selected by further functional and informational constrains. Researchers hope to detect those targets, even they were ultra-trace comparing to major components (source molecules and undefined byproducts). Analytical point of the views, **it is the worst situation.**

Molecular evolution toward Protometabolism¹⁾



1) Aono et al., 2015, Orig. Life Evol. Biosph., 45, 327-338.


Analytical challenges

- Qualification of unpredictable products
70eV EI fragmentation combined with NIST data base was the best way to get structural information => **GC-MS**
- Messy matrix (metal/buffer)
Pretreatment is must before GC-MS:
 - Desalt and metal elimination
 - Molecular screening (chemical prop./ size)
 - Chemical derivatization
 } requiring much time and effort


=> An automated technique using a large volume SPME (**SPME Arrow™**) mounted on a **CTC-PAL system** was applied.

Prebiotic product analyses demand new analytical approaches.

Simultaneous analysis of Potential products in Formose reaction (1861)²⁾



Conventional derivatization procedure³⁾



2) Butlerov, A. C. R., 1861, Acad. Sci., 53, 145-147 3) Borate Minerals Stabilize Ribose A. Ricardo, et al. 2004, Science, 303, 196

SPME Arrow™ extraction and on fiber derivatization

SPME Arrow™ fibers

Conventional SPME	SPME Arrow™
Sorption phase: 0.6 uL	15.3 uL
• Sensitivity	• Sensitivity
• throughput	• throughput
• Robustness	• Robustness
• Reproducibility	• Reproducibility

The simplest schematic workflow

Fiber choice

Polyacrylate (PA)	Carbon WRR
Absorption of polar compounds	Adsorption of diverse molecules having MW<300
Sugars	Amino acids and analogues
	AAA

Standard experiments

GCMS: SHIMADZU QP-2010 Ultra split/splitless injector GC2010 system DB-5 column (Agilent, 30 m, 0.25 mm I.D., 0.25 mm df).


GC conditions		EI-MS conditions		Retention times and reference ions (standard)					
Sugars		AAA		Sugars		AAA			
Injection Temperature	270 °C	Injection Temperature	205 °C	Standard Peaks	Standard sugars	Reference peaks (m/z)	Standard Peaks	Standard amino acids	Reference peaks (m/z)
Split Flow	64.3 ml/min	Split Flow	64.3 ml/min	5.595 2TMS-Glycolaldehyde (C2)		117 116 191 73 59	6.5402TBOMS-pyrvate		73 147 189 261
Injection Mode	Splitless	Injection Mode	Splitless	9.415 3TMS-Glyceraldehyde (C3)		73 147 191 218 307 103	6.84 2TBOMS-Glycine		73 147 218 206
Carrier Gas, Column Flow	He, 1.20 ml/min	Carrier Gas, Column Flow	He, 1.20 ml/min	9.570 3TMS-Glyceraldehyde (C3)		73 103 147 217 307 103	6.760 2TBOMS-Alanine		73 147 158 232 260 302
Oven Temperature	65 °C hold 2 min 1°C/min to 100°C 100°C/min to 280°C	Oven Temperature	65 °C hold 4 min 27°C/min to 280°C 180°C/min to 320 °C hold 2min	8.5854TMS-Ribose (C5)		73 217	7.2202 TBOMS-beta Alanine		73 145 147 218 260 302
		Interface Temperature	250 °C	8.705 4TMS-Ribopyranose (C5)		73 147 191 204 217	8.6403 TBOMS-Serine		73 147 288 302 362 390 432
				8.840 4TMS-Ribofuranose (C5)		73 147 191 204 217			
280°C	200°C/min to hold 2min			8.4304TMS-Arabinose (C5)		73 147 191 204 217			
320°C	200°C/min to hold 2min			8.6704TMS-Arabinopyranose (C5)		73 147 191 204 217			
Interface Temperature	210 °C			8.865 4TMS-Arabinofuranose (C5)		73 147 191 204 217			
				7.430 3TMS-Deoxyribopyranose (C5)		73 75 101 103 116 129 133 147 191 204 217 260			
				7.550 3TMS-Deoxyribose (C5)		73 75 101 103 116 129 133 147 191 204 217 260			

PAL tasks and typical parameter settings

Fiber baking run (background check)		Measurement run		GC-MS data acquisition				
Mount Arrow tool		Fiber exposure in GC inset		Sample extraction	drying	derivatization	desolvation	injection
T [°C]	Sugars	Polyacrylate 100 μm	270	Room temperature	60	45 DMF/HDMS&TMSC=1/1	60	280
	AAA	Carbon WR/PDMS 120 μm	295			60 MTBSFTA/CH3CN=1/1		295
Duration [min]	Sugars		2	60	0.5	10	10	4
	AAA		4	20		20		


Swelling character of PDMS based fiber

The fiber was broken in the derivatization step



derivatization
<heated reagent>

Unexpected swelling of PDMS contained in Carbon WR fiber caused serious damage to the fiber in the derivatization step. It was mechanical problem, so wide sleeve fiber was applied.



Fiber was swelled by the reagents in derivatization.

The expanded surface area would enhance the capacity retaining molecules.⁴⁾ Chemical property of PDMS (log Ko/w (a logarithm of octanol-water partitioning coefficient) > 3.0)⁵⁾ can be modified with appropriate solvents. Wide sleeve fiber is profitable to various applications.


- 4) Ochiai N., Samamoto K. David F., Sandra P. *J Chromatogr A*. 2016, 15, 1455, 45-56.
5) J. Li, Y.-B. Wang, K.-Y. Li, Y.-Q. Cao, S. Wu, L. Wu *Trends Anal. Chem.*, 2015, 72, 141-152.

γ (60Co)-radiolysis sample analyses

The radiolysis sample was prepared by aqueous CH₂O (33 mM) and CH₃CN (30%)/(15%) with NH₄HCO₃. Dose rates of gamma radiation from a ⁶⁰Co source at the Tokyo Institute of Technology's ⁶⁰Co Gamma Radiation Facility, and the control vials were kept in a darkened box at room temperature as a control.

Preliminary product search


direct injection to GCMS Rtx-WAX column



- Similarity search
- 1.2.315 / - / - / - / -
 - 2.2.63 / - / - / - / -
 - 3.3.003 / - / - / - / -
 - 4.3.41 / - / - / - / -
 - 5.3.51 / - / - / - / -
 - 6.3.605 / - / - / - / -
 - 7.3.669 / - / - / - / -
 - 8.4.095 / 4.115 / 4.125
 - 9.4.535 / - / - / - / -
 - 10.4.555 / 4.595 / 4.611
 - 11.5.05 / - / - / - / -
 - 12.5.585 / 5.615 / 5.615
 - 13.6.485 / - / - / - / -
 - 14.6.825 / 6.865 / 6.87
 - 15.7.717 / 7.795 / 7.79
 - 16.8.19 / 8.22 / 8.37
 - 17.8.395 / 8.44 / 8.44
 - 18.10.56 / - / - / - / -
 - 19.12.95 / 12.965 / 12.985
 - 20.13.04 / - / - / - / -
 - 21.13.14 / - / - / - / -
 - 22.13.92 / 13.98 / 13.98
 - 23.14.475 / - / - / - / -
 - 24.14.775 / - / - / - / -
 - 25.16.645 / 16.655 / 16.675
 - 26.19.04 / 19.085 / 19.09
 - 27.19.215 / - / - / - / -
 - 28.19.44 / - / - / - / -
 - 29.19.685 / - / - / - / -
 - 30.20.325 / 20.34 / 20.375
 - 31.20.735 / - / - / - / -
 - 32.21.325 / - / - / - / -
 - 33.21.43 / - / - / - / -
 - 34.22.2 / - / - / - / -
 - 35.15.03


AA analysis using conventional HPLC fluorescence post column derivatization (OPA)

Detected Peaks	Amino acids standard mixture
13.2 uncharacterized	7.7 L-Aspartic Acid
13.6 uncharacterized	6.6 L-Threonine
24.6 uncharacterized	9.4 L-Serine
25.3 uncharacterized	11.1 L-Glutamic Acid
25.8 uncharacterized	13.2 Glycine
27.7 uncharacterized	13.6 L-Alanine
28.1 uncharacterized	14.1 L-Cysteine
30.4 NH ₄ ⁺	15.2 L-Valine
31.0 uncharacterized	17.6 L-Methionine
35.1 uncharacterized	20.5 L-Isoleucine
35.7 uncharacterized	22.4 L-Leucine
	24.3 L-Tyrosine
	25.0 L-Phenylalanine
	28.4 L-Lysine
	29.8 L-Histidine
	30.4 NH ₄ ⁺
	35.1 L-Arginine




Sensitive analyses v.s. Memory effect

Recovery check for 10 min extraction 10 min extraction(60oC)




Signals were substantially decreased in the 2nd analyses for c.a.10 mM samples, indicating that 2 μmol NET AA in 200 μL aqueous solutions can be accumulated for 10 min.




Longer duration in the extraction enhanced signal intensities, and required longer duration in derivatization. Imbalance between long extraction and short derivatization cause serious memory effect.

Comparison between conventional and this technique for 200 uL sugar solution samples.



SPME extraction and on-fiber derivatization



Potential radiolysis products were detected for both 30% acetonitrile and 15% acetonitrile + 1M NH₄HCO₃. Potential TBDMS-AA, Ala and Gly, were detected only in 15% acetonitrile + 1M NH₄HCO₃ sample. Using specific fragment to derivative structures at m/z 73 (C₄H₈Si⁺), we can search the potential target having similar structures to those of amino acids. Carbon WR seemed to eliminate smaller and semi-polar organics detected in direct GC injection using Rtx WAX column. For γ-ray exposure experiment on 33mM formaldehyde, no sugars were detected. The sugar syntheses in radiolysis has not been confirmed. Sugars potentially break under high energy system. This procedure would work in the first screening of the stochastic reaction products. I will apply this technique to not only prebiotic synthetic samples but also biological samples.

Acknowledgements

Authors would like to thank Ms. Shibamoto and Dr. Yoda very much for their support in experiments. Dr. Albert Fahrenbach provided his sample to examine this procedure. This work was partially supported by the Ministry of Education, Culture, Sports, Science and Technology, Japan through a KIBAN(C), JSPS KAKENHI Grant-in-Aid for Scientific Research on Innovative Areas "Hadean Bioscience," Grant Number JP26106003. The radiolysis sample was prepared in the Tokyo Institute of Technology's ⁶⁰Co Gamma Radiation Facility.