

Enzymatic fragmentation vs sonication for the production of cfDNA reference standards for liquid biopsy assays

Aldo Mele¹, Hannah Child¹, Katarzyna Wilczynska¹, Christian Paar² and Julie Wickenden¹
¹ Horizon Discovery, ² Kepler University Hospital, Austria

Introduction

Liquid biopsies hold great promise to revolutionize the field of clinical oncology testing. Cell-free DNA (cfDNA) can be extracted from a routine patient blood sample and used to determine the genetic profile of a solid tumor located elsewhere within the body. This facilitates more informative disease management for the clinician, without the need for invasive surgery for the patient. With new cfDNA NGS assays being able to detect variants from as little as 2-10ng DNA, assay validation to ensure sufficient accuracy has never been so critical. Reference materials that closely mimic real cfDNA samples are essential to support this effort. Here we present results from a comparative study of DNA fragmentation methods applied during the production of cfDNA reference standards. We show a comparison between enzymatic fragmentation and mechanical shearing (sonication), and the benefit of including a size selection step for data accuracy and performance of NGS gene panel workflow.

Methods

DNA extracted from engineered cancer cell lines, representing the Multiplex 1 blend at 5% or 0.1%, was fragmented by mechanical or enzymatic shearing. In addition, a size selection step was included to obtain a fragment size distribution profile that closely mimics real cfDNA samples. The allele frequency of specific variants was confirmed by ddPCR. The eight-sample cfDNA material experimental set was externally tested on the Illumina TruSight Tumor 15 (TST-15) panel to assess library preparation and variant calling performance.

NGS: Sequencing was performed on the MiSeqDx system in RUO mode. Filter settings for analysis with VariantStudio (and automatic analysis) were: Read depth >500 and MAF >2%.

Conclusion


Results show good performance of both sonicated and enzymatically sheared cfDNA material, where all variants present above the LOD could be detected on the Illumina TST-15 assay.

Whilst enzymatically sheared cfDNA did not show any increase in NGS library yield, it did display slightly more accurate variant calling, in addition to a more defined tapestation profile – centred around 168bp, which was further enhanced when combined with a size selection step (Fig 1A yellow trace).

This highlights the potential of these alternative techniques to produce cfDNA that is highly commutable to patient samples and suitable for the validation of ddPCR and NGS liquid biopsy assays.

Results

Tapestation analysis confirmed that both sonication and enzymatic fragmentation produced cfDNA with a peak size in line with real clinical samples (Fig 1). Proof of principle variant detection by ddPCR confirmed the presence of three of the eight mutations across 2 of four genes (EGFR, KRAS, NRAS and PIK3CA) at either 0.1% or 5% variant allele frequency (Fig 2). Library preparation using the Illumina TST-15 gene panel showed good library yield across all eight samples, as assessed by Qubit (Fig 3). NGS sequencing results showed good and comparable variant calling ability between both sonicated and enzymatically sheared samples (Fig 4). (Due to the Limit of Detection (LOD) of the assay used, no variants were detected in any sample at 0.1% AF (even with filter settings removed, data not shown).)



B)

Sample	Region [35-950bp]		
	Peak (bp)	Average Fragment Length (bp)	% of Total
Sonication	178	170	95.08
Sonication + SS	171	180	96.54
Enzymatic	187	243	94.95
Enzymatic + SS	168	196	95.58

Figure 1: Size Distribution
(A) TapeStation analysis of size distribution of all sample conditions (B) Table listing peak size and average fragment length

Library Pool A		ng/μL
1	5% sonic. _ 31097	27.6
2	5% sonic. + SS _ 31168	49.3
3	5% enzymatic _ 31116	44.7
4	5% enzymatic + SS _ 31119	38.5
5	0.1% sonic. _ 29637	50.5
6	0.1% sonic. + SS _ 31169	49.7
7	0.1% enzymatic _ 31117	42.7
8	0.1% enzymatic + SS _ 31120	57.2
Library Pool B		ng/μL
1	5% sonic. _ 31097	49.8
2	5% sonic. + SS _ 31168	53.1
3	5% enzymatic _ 31116	44.5
4	5% enzymatic + SS _ 31119	40.5
5	0.1% sonic. _ 29637	58.7
6	0.1% sonic. + SS _ 31169	52.5
7	0.1% enzymatic _ 31117	45.1
8	0.1% enzymatic + SS _ 31120	47.8

Figure 3: NGS Library Yield

	Mutation	5% AF samples		0.1% AF samples			
		Expected AF (%)	Acceptance Criteria (%)	Measured AF (%)	Expected AF (%)	Acceptance Criteria (%)	Measured AF (%)
Sonication	EGFR L858R	5	3.5 - 6.5	4.5	0.10	0.05-0.15	0.13
	EGFR V769-D770insASV	5	3.5 - 6.5	4.4	0.10	0.07 - 0.20	0.11
	NRAS Q61K	6.3	4.4 - 8.2	6.3	0.13	0.07-0.20	0.14
Sonication + SS	EGFR L858R	5	3.5 - 6.5	4.6	0.10	0.05-0.15	0.13
	EGFR V769-D770insASV	5	3.5 - 6.5	4.2	0.10	0.07 - 0.20	0.10
	NRAS Q61K	6.3	4.4 - 8.2	6.4	0.13	0.07-0.20	0.12
Enzymatic	EGFR L858R	5	3.5 - 6.5	4.7	0.10	0.05-0.15	0.12
	EGFR V769-D770insASV	5	3.5 - 6.5	4.5	0.10	0.07 - 0.20	0.10
	NRAS Q61K	6.3	4.4 - 8.2	6.6	0.13	0.07-0.20	0.13
Enzymatic + SS	EGFR L858R	5	3.5 - 6.5	4.7	0.10	0.05-0.15	0.13
	EGFR V769-D770insASV	5	3.5 - 6.5	4.3	0.10	0.07 - 0.20	0.10
	NRAS Q61K	6.3	4.4 - 8.2	6.5	0.13	0.07-0.20	0.14

Figure 2: Representative ddPCR QC analysis on three of the eight mutations at either 0.1% or 5% variant AF

Expected Variant AF (%) (ddPCR)	Mutations							Sonication (31097)	Sonication + Size Selection (31168)	Enzymatic (31116)	Enzymatic + Size Selection (31119)
	Gene	Variant	Type	Protein Position	AA Change	Exon	Consequence	Variant AF (%) as detected on Illumina TST15 NGS gene panel			
5.0	EGFR	ΔE746 - A750	deletion	0	KELREA/K	19/28	inframe-deletion	7.1	8.2	4.7	5.5
5.0	EGFR	L858R	snv	858	L/R	21/28	missense	4.4	2.7	5.9	6.1
5.0	EGFR	T790M	snv	790	T/M	20/28	missense	5.6	3.6	5.9	4.2
5.0	EGFR	V769 - D770insASV	insertion	766	M/MASV	20/28	inframe-insertion	5.9	4	4.3	3.8
6.3	KRAS	G12D	snv	12	G/D	2/6	missense	7.7	7.9	6.5	5.8
6.3	NRAS	Q61K	snv	61	Q/K	3/7	missense	9.2	5.2	6.5	6.8
6.3	NRAS	A59T	snv	59	A/T	3/7	missense	7	8.3	6.2	6.9
6.3	PIK3CA	E545K	snv	545	E/K	10/21	missense	5.3	6.3	6.5	6.7

Figure 4: NGS results from 5% AF cfDNA test samples run on the Illumina TST-15 assay