


# Applications of chemically modified synthetic guide RNA for CRISPR-Cas9 genome editing

Megan Basila, Eldon Chou, Emily Anderson, Melissa L. Kelley, Anja van Brabant Smith | Dharmacon, part of GE Healthcare, 2650 Crescent Drive, Lafayette, CO 80026, USA

## Abstract

The bacterial CRISPR-Cas9 system has been applied in mammalian cells to efficiently disrupt genes through the formation of targeted DNA double-strand breaks. The Cas9 nuclease is directed to DNA using a guide RNA (gRNA), either as the native dual-RNA system consisting of a DNA-targeting CRISPR RNA (crRNA) and trans-activating crRNA (tracrRNA), or a chimeric single guide RNA (sgRNA) created through the fusion of crRNA and tracrRNA. DNA-free genome engineering can be achieved by using Cas9 mRNA or Cas9 protein with a gRNA, including *in vitro* transcribed (IVT) gRNA, synthetic sgRNA or synthetic crRNA:tracrRNA. While IVT sgRNAs can elicit an immune response, synthetic sgRNA or crRNA:tracrRNA have little to no effect on the immune response and permit chemical modifications to be incorporated to the RNA for increased stability. Here we present chemical modification of synthetic crRNA:tracrRNA with one to three 2'-O-methyl and phosphorothioates (MS) on the 5' and/or 3' ends. These modified RNAs were co-delivered into cells with Cas9 mRNA or Cas9 protein using electroporation. Some modification patterns were found to significantly improve CRISPR-Cas9 gene editing when used with Cas9 mRNA compared to the unmodified versions, yet most modifications did not significantly increase gene editing when used with Cas9 protein. Transfection reagent-mediated delivery of these modified gRNAs into a Cas9-expressing cell line resulted in similar editing efficiencies as the unmodified synthetic gRNAs, and cellular toxicity was observed with certain modification patterns. Of the modifications that were nontoxic, some patterns showed modest improvement in editing efficiency when co-transfected with Cas9 mRNA or Cas9 protein. Overall, our results indicate that MS modifications are required for experiments with co-electroporation of Cas9 mRNA and synthetic gRNA, yet have no impact on editing efficiency when delivered with lipid-based transfection reagents.


## Synthetic guide RNAs show comparable editing to *in vitro* transcribed sgRNA, but no immune response



## Synthetic crRNA:tracrRNA for DNA-free CRISPR-Cas9 gene editing

### Benefits of Dharmacon™ Edit-R™ synthetic guide RNA

- Arrives ready to use (no cloning, sequencing, etc.)
- Transient, fewer off-target effects, less toxic
- Enables high-throughput applications (e.g., arrayed screening)
- Permits application of chemical modifications




## General electroporation workflow for synthetic crRNA:tracrRNA & Cas9 mRNA


### Sequential delivery

### Co-delivery


A protocol for electroporation of Cas9 mRNA and synthetic crRNA:tracrRNA was developed for K-562 cells. One day before electroporation,  $6 \times 10^6$  cells were plated in a 150 mm dish. For sequential electroporations,  $2 \times 10^6$  cells were collected and electroporated with Edit-R Cas9 Nuclease mRNA (Cat #CAS11195; 5  $\mu$ g) using the Lonza Nucleofector 2b™, as per the manufacturer's protocol. Electroporated cells were plated for 6 hours, collected and electroporated with crRNA:tracrRNA (5.4  $\mu$ M). Cells were plated again and incubated for 72 hours and analyzed for gene editing. For co-delivery, K-562 cells were electroporated with Cas9 mRNA and crRNA:tracrRNA, as described above, in a single electroporation. Electroporated cells were plated and incubated for 72 hours then analyzed for gene editing.




## Sequential electroporation is required for delivery of unmodified crRNA:tracrRNA & Cas9 mRNA




## Modifications of synthetic guide RNAs for increased stability




## Modification of both ends of crRNA:tracrRNA stabilize RNAs for co-electroporation with Cas9 mRNA




## Modification of the single-stranded regions of crRNA:tracrRNA is sufficient for stabilization in co-electroporation with Cas9 mRNA




## Modification of crRNA:tracrRNA can increase gene editing in co-electroporation with Cas9 protein



## RNA modifications do not affect gene editing activity with lipid transfection in Cas9-expressing cells, but some patterns negatively affect cell viability



## Some modification patterns show modest improvements in gene editing with co-transfection of modified crRNA:tracrRNA & Cas9 mRNA or Cas9 protein



## Conclusions

- Stabilizing modifications on crRNA:tracrRNA are required for co-electroporation with Cas9 mRNA
- Some modifications improve gene editing efficiency in co-electroporation with Cas9 protein for some gene targets
- Lipid transfection of crRNA:tracrRNA in a Cas9 stable cell line shows no difference in gene editing between unmodified and modified RNAs
- Some modification patterns are toxic to cells with lipid transfection
- Stabilization of the single-stranded regions of the dual RNAs modestly increases gene editing in lipid co-transfection with Cas9 mRNA or Cas9 protein

[gelifesciences.com/dharmacon](http://gelifesciences.com/dharmacon)

GE, the GE monogram and GE Healthcare are trademarks of General Electric Company. Dharmacon, Inc. is a General Electric company doing business as GE Healthcare. All other trademarks are the property of General Electric Company or one of its subsidiaries. ©2016 General Electric Company. All rights reserved. First version published December 2016. GE Healthcare UK Limited, Amersham Place, Little Chalfont, Buckinghamshire, HP7 9NA, UK