




# Interfacing Microwave Synthesis with Enabling Technologies in Drug Discovery



Doris Dallinger and C. Oliver Kappe\*

Institute of Chemistry, Organic and Bioorganic Chemistry, Karl-Franzens University Graz, Heinrichstrasse 28, A-8010-Graz, Austria  
email: do.dallinger@uni-graz.at website: http://www.maos.net

## Introduction

Modern drug discovery steadily relies on high-speed organic synthesis and combinatorial chemistry techniques for the rapid generation of compound libraries. Microwave-assisted organic synthesis in combination with polymer-assisted solution phase (PASP) methods are powerful techniques that are frequently applied in both academia and industry.


Here we report on the rapid synthesis/decoration of privileged heterocyclic scaffolds, namely the dihydropyrimidine (DHPM) scaffold, by applying microwave heating in conjunction with PASP approaches. The use of resin-bound catalysts, reagents and scavengers as well as fluororous reagents in conjunction with microwave heating will be described for amidations, esterifications *via* a Mitsunobu protocol [1] and *N*-acylations [2].

[1] Desai, B.; Dallinger, D.; Kappe, C. O. *Tetrahedron* **2006**, 62, in press.

[2] Dallinger, D.; Gorobets, N. Yu.; Kappe, C. O. *Org. Lett.* **2003**, 5, 1205. Dallinger, D.; Gorobets, N. Yu.; Kappe, C. O. *Mol. Diversity* **2003**, 7, 229.


## 1 Polymer-supported catalysts

Acylation is used for a protection purpose at the *N*3-position of the DHPM. For ease of purification, polymer-bound DMAP is used as acylation catalyst.



## 2 Polymer-supported scavengers

In order to synthesize a library of *N*3-acylated DHPMs, we employed this high-throughput protocol which includes a scavenging step to remove both excess anhydride and bis-acylated byproduct.





➤ Several scavenging reagents **A** were evaluated both under rt and MW (80-100 °C) conditions (see data **B** and **C**):

- A**
- Polystyrene-bound ethylenediamine
  - Functionalized silica gel ethylenediamine
  - StratoSpheres Plugs (diethylenetriaminomethyl)
  - SynPhase Lanterns (aminomethyl)




Room Temperature (25 °C)

Microwave Heating (80-100 °C)




## 3 Polymer-supported reagents

Polymer-bound carbodiimide with subsequent SPE purification is used in this amidation protocol at the DHPM-C-5 position.



## 4 Fluororous reagents

To overcome the laborious purification of the Mitsunobu reaction, we decided on a fluororous approach using fluororous Mitsunobu reagents for the esterification at the DHPM-C-5 position.



## 5 Conclusion

The combination of MW heating with PASP methods is a powerful tool for organic synthesis. Not only an acceleration of reactions is possible but also the purification issue is simplified. Here we showed examples for the use of polymer-bound catalysts, scavengers and reagents as well as fluororous reagents to ease the sometimes laborious reaction work-up.

## Acknowledgements

This work was supported by the Austrian Science Fund (FWF, I18-N07). We thank Biotage (Uppsala, Sweden) for the use of the Emrys Synthesizer and Initiator Eight. We also thank Silicycle Inc., Mimotopes and Polymer Laboratories Ltd. for providing the scavengers used in this study and FTI for the fluororous Mitsunobu reagents.