

# Hepatitis B virus (HBV) and Human immunodeficiency virus (HIV) antibodies detected by peptide microarrays

Ahmed Abd El Wahed<sup>1</sup>, Ulrike Beutling<sup>2</sup>, Ronald Frank<sup>2</sup>, Gerhard Hunsmann<sup>1</sup>, and Hans-Joachim Fritz<sup>3</sup>

1) Institute of Virology, Universitätsmedizin Göttingen, Georg August University Göttingen,

2) Department of Chemical Biology, Helmholtz Centre for Infection Research, Braunschweig,

3) Department of Molecular Genetics and Preparative Molecular Biology, Institute of Microbiology and Genetics, Georg August University Göttingen.

## Introduction

- ✓ The identification of HBV and HIV epitopes is important for the development of novel diagnostics and vaccines
- ✓ In this study, we have developed HBV and HIVenv chips with overlapping oligopeptides encompassing the full amino acid sequences of different HBV and HIV polypeptides. In addition, a random peptide library composed of 4608 15-mers was prepared.
- ✓ The chips were used for analyzing monoclonal antibodies and sera from HIV and HBV infected individuals.

## Results

### 1. Targets of monoclonal antibodies




Figure 1: Target sequences of four monoclonal antibodies to HIVgp160 comprise 12, 9, 6, or 3 amino acids.

### 2. Analysis of sera from patients recovering from HBV infection and of human sera neutralizing HIV-1

Sequence recognized	Position
NSNNPD	aa 37-43 of preS1, genotype A
LGFFPD PLGFFP	aa 22-27 of preS1, genotype A aa 10-16 of preS1, genotype D
PHGGVLGWSPQA PPHGGLLGW	aa 69-84 of preS1, genotype A aa 58-66 of preS1, genotype D
DPKVRGLYF RVRGLY	aa 13-21 of preS2, genotype A aa 15-20 of preS2, genotype D
PISSIFSRIGD	aa 39-49 of preS2, genotype D
PYKMDIDPY	aa 7-15 of HBcAg, genotype A
SVRDLLDNASALYRE PSVRDLDL TASALYR	aa 26-40 of HBcAg, genotype D aa 25-39 of HBcAg, genotype A

Table 1: Epitopes detected by 5 sera of patients recovering from HBV infection.



Figure 2: Three epitopes detected with two HIV-1 neutralizing human sera.

### 3. Screening the random peptide library with monoclonal antibody MA18/7




Figure 3: MA18/7 detected 10 strong responders (white squares) and 7 weak responders.

Q31B08	L	Q	L	Q	G	Y	M	H	R	V	E	W	C	Y	D (7)	P (3)	A (5)	F (8)	
Q38G10	T	Y	W	D	R	G	F	Q	G	W	Y	G	M	I	N	E (3)	A (2)	G (3)	Y
Q17D08	V	D	Q	A	F	K	T	F	Q	A	E	R	K	H	T				
Q13D02	G	R	V	A	L	Y	T	E	P	G	F	R	V	Q	Y				
Q45D11	V	M	P	Y	P	E	T	E	P	A	F	L	D	K	C				
Q31B07	C	W	F	N	C	M	K	M	D	P	G	F	K	T	I				
Q24G04	T	M	Q	H	D	K	C	W	Q	Y	W	F	L	C	S				
Q45T01	S	A	Y	Q	I	F	M	E	M	W	S	D	K	A	F				
Q43C04	M	A	S	E	F	T	Q	A	L	D	A	A	F	F	K				
Q44H09	S	N	F	P	D	A	A	F	N	N	Q	E	G	I	D				

Consensus: DPAF

Table 2: Sequence of the ten strongly responding peptides. Common motif highlighted in yellow.

Table 3: The profile



Figure 4: Verification of the target sequence of mAb MA18/7 using HBV scanning chip.




Figure 5: Single residue substitution analysis revealing the range of residues permissible at each motif position.

## Conclusions

- ✓ HBV and HIVenv chips are powerful tools to identify and map humoral immune responses against HBV and HIV.
- ✓ The random peptide library has proven potential for identifying B-cell epitopes without prior knowledge of immunizing antigen.
- ✓ Using the random peptide library, studies have been initiated to identify discontinuous epitopes.

Acknowledgments:  
We thank Dr. Angela Uy for HBV sera and Programme EVA (NIBSC, UK) for HIV-1 mAbs and sera neutralizing HIV-1 in cell culture