




MICROWAVE-ASSISTED SCAVENGING OF ELECTROPHILES UTILIZING POLYMER-SUPPORTED SEQUESTRATION REAGENTS



Doris Dallinger and C. Oliver Kappe*

Institute of Chemistry, Organic and Bioorganic Chemistry, Karl-Franzens University Graz, Heinrichstrasse 28, A-8010-Graz, Austria
email: do.dallinger@uni-graz.at
website: http://www.maos.net


1 Introduction

In recent work we have reported the automated generation of a diverse library of multifunctionalized dihydropyrimidines (DHPMs) utilizing a microwave-mediated solution phase Biginelli three component condensation [1]. Since most of the pharmacologically attractive DHPM derivatives are *N3*-acylated analogs [2], we became interested in developing a rapid method for accessing libraries containing this structural motif in high-throughput format. We now present the high-speed scaffold decoration of this library, introducing a new point of diversity at the *N3* position [3]. Different scavenging techniques using polymer-supported sequestration agents are described for the purification steps in the synthesis of *N3*-acylated dihydropyrimidines.

- [1] Stadler, A.; Kappe, C. O. J. *Comb. Chem.* **2001**, *3*, 624
[2] Kappe, C. O. *Acc. Chem. Res.* **2000**, *33*, 879
[3] Dallinger, D.; Gorobets, N. Yu.; Kappe, C. O. *Org. Lett.* **2003**, *5*, 1205

2 Microwave-Assisted *N3*-Acylation

Selective Acylations using Scavenging Techniques



Dallinger, D.; Gorobets, N. Yu. Kappe, C. O. *Mol. Diversity* **2003**, *7*, 229

3 Kinetic Analysis of Excess Bz_2O Scavenging using Supported Sequestration Reagents

In order to make this protocol amenable to a high-throughput format, several purification issues needed to be considered:


- Excess anhydride
- Acid and DMAP
- Bis-acylated Byproduct

Several scavenging reagents **A** were evaluated both under rt and MW (80-100 °C) conditions (see data **B** and **C** on the right):


- Polystyrene-bound ethylenediamine
- Functionalized silica gel ethylenediamine
- StratoSpheres Plugs (diethylenetriaminomethyl)
- SynPhase Lanterns (aminomethyl)

Complete sequestration of Bz_2O was also achievable (MW: 10 min/100 °C) with Lanterns by mimicking the reaction conditions (without DHPM) but at a higher level of scavenger concentration (2.4 equiv.).


Automation was possible using water as scavenger at 100 °C (removing excess anhydride and bis-acylated product) coupled with an SPE purification protocol (removing acid and DMAP, see HPLC chromatograms at 280 nm **D** on the right)



Room Temperature (25 °C)



Microwave Heating (80-100 °C)



4 Conclusion

- Kinetic investigations of 4 different polymer supported amine scavengers
- Reducing scavenging time from hours to minutes using microwave heating
- 44 member library of acylated DHPMs (80 % average yield)
- Shorter reaction times utilizing a microwave protocol
- Liquid handling possible using water as sequestration agent coupled with SPE

Acknowledgements

This work was supported by the Austrian Science Fund (FWF, P15582). We thank PersonalChemistry AB (Uppsala, Sweden) for the use of their Emrys Synthesizer. We also thank Silicycle Inc., Mimotopes and Polymer Laboratories Ltd. for providing the scavengers used in this study.