

DC microelectrode array for cell investigation

Marisa Rio¹, Sharanya Bola², Salvatore Girardo³, Richard H. Funk², Gerald Gerlach¹

¹ Institute of Solid State Electronics, Faculty of Electrical and Computer Engineering, Technische Universität Dresden, Germany

² Department of Anatomy, University Hospital Carl Gustav Carus, Technische Universität Dresden, Germany


³ Biotechnology Center, Technische Universität Dresden, Germany

Abstract


- Endogenous DC electric fields are known play a role in many cell biological phenomena, ranging from cell adhesion and migration, embryonic and tissue development to wound healing;
- Recent studies suggest that electrically sensitive cells lines are the rule and not the exception. However, the lack of understanding of these phenomena conditions mimicking and optimization;
- In this work, a microelectrode array enabling the use of DC current was designed. The first results shown the device is biocompatible;
- Future work will systematically investigate the effect of DC electric fields on 2D cultures of 661w murine photoreceptor-derived cells.

Microelectrode array fabrication


Step 1: Photolithographically patterned SU-8 master




Step 2: PDMS replica by soft lithography using the MicCell platform (GeSiM GmbH)




Step 3: Photolithographically patterned SU-8 free-standing layer



Step 4: Surface modification



Step 5: Multilayer integration



Results

The device's biocompatibility was first tested by accessing cell growth and adherence of the 661w cells on the (3-Aminopropyl)triethoxysilane (APTES) treated SU-8 membranes.

Retinal cell culture of murine 661W photoreceptor-derived cell were seeded with 8×10^4 concentration and allowed to grow for 48 hours. In order to access cell number and proliferation in the different as-treated and protein coated membranes, the DNA-intercalating fluorochromes, Hoechst 33342 and Hoechst 33258 were used. Cells were subsequently washed with PBS, dried, mounted with DABKO solution and immediately imaged.




Figure 1: 661w stained with Hoechst growing on a) APTES-treated SU-8 membrane; b) APTES-treated SU-8 membrane + Fibronectin

The cultures of 661w cell showed to adhere and grow on the APTES-treated SU-8 membranes independently of the addition of the glycoprotein fibronectin.

Conclusions and Future Work

- In order to mimic the endogenous electric fields in vivo the presented DC MEA uses microfluidic channels to separate the electrodes from the cell surface
- The DC MEA uses multilayer integration to construct a platform that enables application of DC electric fields for cell investigation under an inverted microscope
- 661w cells showed good adhesion and growth on the APTES-treated SU-8 membranes
- The use of fibronectin to improve cell adhesion was showed unnecessary for APTES-treated SU-8 membranes
- DC electric field stimulation and recording under different conditions including addition of growth factors will be the next step

Bibliography

- [1] Du Bois-Reymond, E. Vorläufiger Abriss einer Untersuchung über den sogenannten Froschstrom und die electromotorischen Fische. 1843, Ann. Phy. U. Chem. 58, pp. 1-30.
- [2] Galvani L. Aldini, G. Aloysii Galvani: de viribus electricitatis in motu musculari. 1792, Apud Societatem Typographicam.
- [3] McCaig, C. Rajnicek, A.M. Song, B. Zhao, M. Controlling cell behavior electrically: current views and future potential. 2005, Physiol. Rev. 85, pp. 943-978.
- [4] Funk, R. H.W. Monsees, T. Özkucur, N. Electromagnetic effects – From cell biology to medicine. 2009, Progress in histological and cytochemistry 43, pp. 177-264.
- [5] Messerli, M. A. Graham, D.A. Extracellular Electrical Fields Direct Wound Healing and Regeneration. 2011, Biol. Bull. 221, pp. 79-92.
- [6] Levin, M. Large-scale biophysics: ion flows and regeneration. 2007, Trends in Cell Biology, pp. 261-270.
- [7] Derix, Jonathan. Biokompatibles Gleichstrom-Mikroelektrodenarray. Dresden : TUD Press, 2010.

Acknowledgments

This work is part of the research program of the Research Training group "Nano and Bio Technologies for the Packaging of Electrical Devices" financially supported by the German Research Foundation (DFG, GRK 1401).